**H&D srl**Strada Langhirano 264/1A - Parma (PR) - ITALY  
+39 0521 462607 - info@hedsrl.it

## SAT TEST

**REF** SAT50

50 test

**REF** SAT100

100 test

### PRINCIPLE OF THE METHOD

The defense against damaging attack at the mouth level caused by reactive species, in particular by free radicals, is guaranteed by the saliva, which contains the main constituents of the antioxidant barrier. It includes endogenous and exogenous substances that have their own "antioxidant capacity", in relation to their redox potential. The basic principle of the SAT test exploits the ability of a colored solution of ferric ions (Fe<sup>3+</sup>) complexed to a chromogen, to bleach, when Fe<sup>3+</sup> ions are reduced to ferrous ions (Fe<sup>2+</sup>) in the presence of reducing agents (antioxidants).

### KIT COMPONENTS

	SAT50	SAT100
<b>Reagents:</b>		
<b>- Reagent R1 SAT TEST</b>		
Chromogenic mixture in the cuvette, pre-dosed. Store at 15-25°C	50 pcs	100 pcs
<b>- Reagent R2 SAT TEST</b>		
Ferric nitrate solution, stabilizers and preservatives. Store at 15-25°C	1x3ml	2x3 ml
<b>Materials:</b>		
- Dosed cuvette of reagent R1, disposable/single-dose	50 pcs	100 pcs
- vial of reagent R2;	1 pcs	2 pcs
- tips, disposable/single-dose.	100 pcs	200 pcs
- 50 cotton squares	50 pcs	100 pcs
-50 plastic glasses	50 pcs	100 pcs
-50 cards values	50 pcs	100 pcs
- 3 plastic tweezers	3 pcs	6 pcs

### MATERIAL REQUIRED BUT NOT SUPPLIED

Pipette (10 µl and 40 µl) - Photometer

### PRECAUTIONS AND WARNINGS

According to Reg 1272/2008 CLP, the product is not classified as dangerous.

Handle the product with caution, according to good laboratory practice, avoid ingestion, avoid contact with skin, eyes and mucous membranes. Safety data sheets of individual components are available upon request.

Treat all the samples as if they contain infectious agents. Observe precautionary practices against microbiological hazards by following the procedures and standards for proper sample disposal.

Do not use after the expiration date.

### SPECIMEN

Biological fluid: saliva

Storage of the saliva:

0-4°C	-20°C	Room temperature
24 hours	7 days	2 hours

The first test should be performed on patients fasting from the evening before. This is because food and drinks, such as coffee, tea, fruit juice, affect the level of antioxidants in plasma and thus distort the test result. There are no evidences of substantial fluctuations in value over time if the subject conditions remain unchanged. If the patient is not fasting is advisable to have the patient drink a glass of water (100-150 mL) half an hour before the test.

### SAMPLE PREPARATION

-Give the piece of cotton square to the patient and invite him to roll it in the mouth for 1 minute (for about 60 bites without chewing the cotton) in order to induce the production of saliva. Saliva should be conveyed on the cotton  
- Squeeze cotton in the glass using the tweezers supplied, eventually help with the wall of the glass.  
-Dispose of the cotton according to applicable regulations and wash the tweezers with water and make sure to dry it with paper towels.

### TEST PROCEDURE

**ATTENTION:** Before carrying out the tests, you should prepare all necessary materials. The instrument must be switched on at least 10 minutes before a test and should not be used until it has finished warming up.

-Take the cuvette containing the SAT R1 reagent and add 40 µl of R2 reagent solution using the pipette equipped with the SAT box and the relative disposable tip.  
-Close the cuvette with the lid and mix by inversion for about 10 seconds  
-Insert the cuvette into the reading cell, making sure that the ribbed sides are oriented according to the instruction reported and in accordance with the label and also make sure that the cuvette is pushed until the bottom of the reading cell. The instrument carries out the first reading in about 2 seconds. Remove the cuvette from the reading cell following the instructions that appear on the display.  
-Take a sample of 10 µl of saliva and add it into the cuvette containing the R1+R2 SAT solution which you just removed from the reading cell of the instrument. The saliva must be taken using the white pipette equipped with the SAT box and the relative disposable tips.  
-Close the cuvette with the lid and mix by inversion for about 10 seconds. Insert the cuvette into the reading cell, as above, making sure that the ribbed sides are oriented according to the instruction reported and in accordance with the label and also make sure that the cuvette is pushed until the bottom of the reading cell. The instrument carries out the second reading. Wait 1 minute for the result and the printing of the receipt with the value (µmol/l)

### REFERENCE VALUES

The reference values are expressed in micromol per liter (µmol/L or µM) of antioxidants. In this case the reference is the antioxidant vitamin C which is used as iron-reducing agent of interest:

< 1000 µM	Severe shortage
1000 – 1500 µM	Optimal normal values
1500 – 2000 µM	Normal values
2000 – 2500 µM	Border line values
> 2500 µM:	Possible inflammatory processes on place

### ANALYTICAL PERFORMANCE

Unit of measure: 1 U Cor = 1.4 µmol / liter of vitamin C

Linearity: the method is linear in the 300-7000 U Cor range

Precision: CV% &lt;5.5%

Interferences: no interference was observed in the presence of a phosphate concentration of less than 40 mg/dl;

### WASTE DISPOSAL

The product should be disposed of according to local legislation on waste management

### BIBLIOGRAPHY

Cornelli U et al. Pan Minerva Medica 2010; 52(2).  
Benedetti S et al. Clin Lab Med 2014;60 475-482

### SYMBOLS

Check the instructions for use

Temperature range in which to store the product

Keep away from direct light sources

Manufacturer

Disposable

Catalog number

Lot number

For in vitro diagnostic use only

Expiry date (year / month)

**Please read this leaflet carefully before proceeding to use the product. No liability is accepted for damage resulting from improper use or use not covered in this package insert.**

### CAUTION: CHOKING HAZARD

**Do not perform the test on subjects under 36 months of age or on subjects who in the opinion of the user, due to age, deficiencies or diseases, are at risk of suffocation during saliva collection**

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