**H&D srl**Strada Langhirano 264/1A - Parma (PR) - ITALY
+39 0521 462607 - info@hedsrl.it

d-ROMs test FAST

REF IKIT50FAST

50 test

REF IKIT100FAST

100 test

PRINCIPLE OF THE METHOD

Reactive oxygen metabolites, referred to as ROMs in scientific literature, are a variety of free radicals which are characterized by the odd number of electrons found around the external orbit of oxygen. Given their extreme chemical instability, they form derivatives, in both plasma and cells, which remain highly chemically reactive and have an effective oxidative capacity. These derivatives, which react with a particular chromogen which has been correctly buffered, form a colored compound which may be measured photometrically with a maximum absorbency peak at 505 nm and directly proportional to their concentration.

KIT COMPONENTS

Reagents:	IKIT50FAST	IKIT100FAST
- Reagent R1 d-ROMs Test FAST Chromogenic mixture condensed in the cuvette, pre-dosed. Store at 15-25 °C	50 pcs	100 pcs
- Reagent R1 d-ROMs Test FAST Chromogenic mixture condensed in the cuvette, pre-dosed. Store at 15-25 °C	50 pcs	100 pcs
- Reagent R3 d-ROMs Test FAST Catalyst solution. Store at 15-25 °C	1x2ml	2x2 ml
Materials:		
- tips, disposable	100 pcs	200 pcs
- sterile lancets, disposable	50 pcs	100 pcs
- heparinized microvettes, disposable	50 pcs	100 pcs

MATERIAL REQUIRED BUT NOT SUPPLIED

Pipette (10 µl) - Photometer

PRECAUTIONS AND WARNINGS

According to Reg 1272/2008 CLP the product is not classified as dangerous.

Handle the product with caution, according to good laboratory practice, avoid ingestion, avoid contact with skin, eyes and mucous membranes. Safety data sheets of individual components are available upon request.

Treat all the samples as if they contain infectious agents. Observe precautionary practices against microbiological hazards by following the procedures and standards for proper sample disposal.

Do not use after the expiration date.

SPECIMEN

Biological fluid: heparinized capillary blood from which we obtain plasma by centrifugation. **Do not use blood treated with citrate or EDTA** or other iron chelators. Only use ethyl alcohol as a disinfectant.

Storage of the sample:

18-25°C	0-4°C	-20°C
24 hours	48 hours	1 year

The test can be performed at any time during the day.

Sensitive variation of the value do not happen during the time if the conditions of the patient do not vary.

SAMPLE PREPARATION

Prepare a microvette by pulling it out from its container and pulling the small cap attached to the main cap.

- Gently massage the fingertip and disinfect it with alcohol (strictly avoid hydrogen peroxide, a powerful oxidant).
- Practice a puncture on a fingertip with a sterile lancet and gently massage the finger in order to promote blood spill.
- Wipe away the first drop of blood (rich in cellular liquid) with a cotton ball.
- Move the finger close to the microvette and let the blood enter through the smaller hole. Fill the microvette until the end of the fins.
- Close the small hole with the cap FIRST and THEN the main cap. Finally insert the microvette into the container.

TEST PROCEDURE

Warning: Before the exams, you should prepare all the necessary material. The instrument must be switched on at least 10 minutes before the start of a test and should not be used until it has completed the warm-up phase.

- Place the microvette in the centrifuge, with the appropriate counterweight, and start the exam. The sample is centrifuged for 90 seconds in order to separate the plasma.
- Prepare the working solution by depositing the micro test tube containing the reagent R2, 10 µl of reagent R3 using the white pipette date supplied with the instrument and shake for reversing for about 10 seconds.
- After centrifugation and preparation of the working solution taken with the white pipette that comes with the tool, 10 µL of plasma and deposited in micro test tube containing the working solution and mix by inversion for at least 10 seconds.
- Transfer the contents of micro test tube, namely the working solution in which the sample was diluted, into the cuvette containing the predosed R1 reagent.
- Close the vial with the cap and mix by inversion for at least 10 seconds.
- Insert the cuvette in the reading room of the tool making sure that the knurled sides are oriented as indicated by the label placed on the label tool. The instrument will perform the analysis in 150 seconds.

REFERENCE VALUES

>500 U Carr	Very high oxidative stress
400-500 U Carr	High oxidative stress
340-400 U Carr	Oxidative stress
320-340 U Carr	Slight oxidative stress
300-320 U Carr	Borderline range
250-300 U Carr	Normal values

ANALYTICAL PERFORMANCE

Unit of measure: 1 U Carr = 0.08 mg / dl hydrogen peroxide

Linearity: the method is linear in the 50-600 U Carr range

Precision: CV% <5.0%

Interference: the addition of anticoagulants capable of chelating iron, such as EDTA or citrates, give rise to underestimates of the data; the use of disinfectants different by ethyl alcohol can lead to anomalous results.

WASTE DISPOSAL

The product should be disposed of according to local legislation on waste management.

BIBLIOGRAPHY

- Cornelli U, et al. Intern. Union of Angiology's Bulletin. 1999. 15: 7-10.
 Cesaroni MR, et al. International Angiology. 1999. 18 (2): 127-130.
 Alberti A, et al. Res Chem Intermed. 2000. 26 (3): 253-67.
 Trotti R., et al. 2001. Haematologica. 86: 85-91.
 Gerardi GM, et al. Clin Chem Lab Med. 2002. 40 (2): 104-110.
 Cornelli U. et al. JCDSA 2011; 1:64-70

SYMBOLS

Check the instructions for use

Temperature range in which to store the product

Keep away from direct light sources

Manufacturer

Disposable

Catalog number

Lot number

For in vitro diagnostic use only

Expiry date (year / month)

Please read this leaflet carefully before proceeding to use the product. No liability is accepted for damage resulting from improper use or use not covered in this package insert.

**H&D srl**Strada Langhirano 264/1A - Parma (PR) - ITALY
+39 0521 462607 - info@hedsrl.it

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