**H&D srl**Strada Langhirano 264/1A - Parma (PR) - ITALY
+39 0521 462607 - info@hedsrl.it

PAT TEST

REF PAT50

50 test

REF PAT100

100 test

PRINCIPLE OF THE METHOD

At the blood level, defence against the lesion caused by the reactive species, and in particular the ROS (reactive oxygen species), is guaranteed by the plasma antioxidant barrier; They are both endogenous and exogenous substances that have their own antioxidant capacity in relation to their oxidation potential. This antioxidant capacity is related to the properties of the barrier components to yield equivalent reducing reagents to the reactive species. The PAT test is based on its ability to have a coloured solution of ferric ions (Fe^{3+}) complexed by a chromogen to decolour when ferric ions are reduced to ferrous ions in the presence of reducing substances. The PAT Test measures the real plasma antioxidant power because excludes potential interfering reactions such as phosphates.

KIT COMPONENTS

	<u>PAT50</u>	<u>PAT100</u>
Reagents:		
- Reagent R1 PAT TEST Chromogenic mixture in the cuvette, pre-dosed. Store at 15-25°C	50 pcs	100 pcs
- Reagent R2 PAT TEST Ferric nitrate solution, stabilizers and preservatives. Store at 15-25°C	1x3ml	2x3 ml
Materials:		
- tips, disposable	100 pcs	200 pcs
- sterile lancets, disposable	50 pcs	100 pcs
- heparinized microvettes, disposable	50 pcs	100 pcs

MATERIAL REQUIRED BUT NOT SUPPLIED

Pipette (10 µl and 40 µl) - Photometer

PRECAUTIONS AND WARNINGS

According to Reg 1272/2008 CLP the product is not classified as dangerous.

Handle the product with caution, according to good laboratory practice, avoid ingestion, avoid contact with skin, eyes and mucous membranes. Safety data sheets of individual components are available upon request.

Treat all the samples as if they contain infectious agents. Observe precautionary practices against microbiological hazards by following the procedures and standards for proper sample disposal.

Do not use after the expiration date.

SPECIMEN

Biological fluid: heparinized capillary blood from which we obtain plasma by centrifugation. **Do not use blood treated with citrate or EDTA** or other iron chelators. Only use ethyl alcohol as a disinfectant.

Storage of the sample:

18-25°C	0-4°C	-20°C
12 hours	48 hours	2 months

The test must be performed on the patient fasting from the previous evening. The food and drinks going to change the level of antioxidant in systematically and therefore distort the test results. Not occur substantial change in the value over if you do not vary the condition of the subject.

SAMPLE PREPARATION

Prepare a microvette by pulling it out from its container and pulling the small cap attached to the main cap.

- Gently massage the fingertip and disinfect it with alcohol (strictly avoid hydrogen peroxide, a powerful oxidant).
- Practice a puncture on a fingertip with a sterile lancet and gently massage the finger in order to promote blood spill.
- Wipe away the first drop of blood (rich in cellular liquid) with a cotton ball.
- Move the finger close to the microvette and let the blood enter through the smaller hole. Fill the microvette until the end of the fins.
- Close the small hole with the cap FIRST and THEN the main cap. Finally insert the microvette into the container.

TEST PROCEDURE

Warning: Before the tests, you should prepare all the necessary material. The instrument must be switched on at least 10 minutes before the start of a test and should not be used until it has completed the warm-up phase.

- Place the microvette in the centrifuge with the appropriate counterweight, and start the test. Centrifuge the sample for 90 seconds in order to separate the plasma.
- Take the cuvette containing the reagent R1 and add 40 µl of reagent R2 by using the green pipette supplied with the instrument and its disposable tips.
- Close the cuvette with the cap and mix by inversion for exactly 10 seconds.
- Insert the cuvette in the reading cell of the instrument making sure that the knurled sides are oriented as indicated by the label on the instrument. The instrument completes the first reading in about 2 seconds. At the end remove the cuvette from the reading cell.
- Add 10 µl of plasma to the solution R1 + R2 contained in the cuvette. Plasma should be taken using the white pipette supplied with the instrument and its disposable tips.
- Close the cuvette and mix by inversion for at least 10 seconds. Insert the cuvette in the reading cell. The instrument completes the second reading in 1 minute.

REFERENCE VALUES

>2800 U Cor	Very high value
2800-2200 U Cor	Normal value
2200-2000 U Cor	Borderline range
2000-1800 U Cor	Slight deficiency status
<1800 U Cor	Deficiency status

ANALYTICAL PERFORMANCE

Unit of measure: 1 U Cor = 1.4 µmol / liter of vitamin C

Linearity: the method is linear in the 500-10,000 U Cor range

Precision: CV% <5.5%

Interferences: no interference was observed in the presence of a phosphate concentration of less than 40 mg/dl; the addition of anticoagulants capable of chelating iron, such as EDTA or citrates, give rise to overestimates of the data; the use of disinfectants different by ethyl alcohol can lead to anomalous results.

WASTE DISPOSAL

The product should be disposed of according to local legislation on waste management

BIBLIOGRAPHY

Cornelli U, et al. *Intern. Union of Angiology's Bulletin*. 1999. 15: 7-10.
 Cesarone MR, et al. *International Angiology*. 1999. 18 (2): 127-130.
 Alberti A, et al. *Res Chem Intermed*. 2000. 26 (3): 253-67.
 Trotti R, et al. 2001. *Haematologica*. 86: 85-91.
 Gerardi GM, et al. *Clin Chem Lab Med*. 2002. 40 (2): 104-110.

SYMBOLS

Check the instructions for use

Temperature range in which to store the product

Keep away from direct light sources

Manufacturer

Disposable

Catalog number

Lot number

For in vitro diagnostic use only

Expiry date (year / month)

Please read this leaflet carefully before proceeding to use the product. No liability is accepted for damage resulting from improper use or use not covered in this package insert.

**H&D srl**Strada Langhirano 264/1A - Parma (PR) - ITALY
+39 0521 462607 - info@hedsrl.it

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