H&D srl

Strada Langhirano 264/1A - Parma (PR) - ITALY
+39 0521 462607 - info@hedsrl.it**REDOX FAST**

d-ROMs Test FAST + PAT Test

REF REDOXBFP 50 test**REF** REDOXBFP20 20 test**PRINCIPLE OF THE METHOD**

In every living organism exists a delicate balance between the production of "free radicals" and their elimination at the hands of antioxidant systems. We call oxidative stress (OS) the breakdown of this balance. A situation of OS derives from an increased production of free radicals or by a reduced responsiveness by of the antioxidant network or in the presence of both phenomena. Free radicals are molecules with an odd number of unpaired electrons in the outermost orbital and generate during specific chemical, biochemical or physical processes starting mainly from reactive oxygen species (ROS). At systemic level, the defense against the detrimental attack by free radicals is guaranteed by the antioxidant network in which there are exogenous substances (eg. vitamins C and E) and endogenous ones (eg. enzymes such as catalase). From an analytical point of view, we may measure the main ROS and the antioxidant defenses with two tests: d-ROMs Test FAST and PAT Test. The d-ROMs Test FAST quantifies the concentration of peroxides (the main ROS in the plasma) by a colorimetric reaction: such analytes placed in a suitable buffer react with the appropriate chromogenic coloring the solution. The coloring of the solution will be proportional to the concentration of the analytes. The PAT Test measures the iron-reducing power by a colorimetric reaction: it is based on the ability to discolor of a colored solution of ferric ions (Fe³⁺) complexed to a chromogen, when the ferric ions are reduced to ferrous ions (Fe²⁺) in the presence of reducing substances. The decolorization of the solution will be proportional to the concentration of antioxidants. The PAT Test measures the real plasma antioxidant power because excludes potential interfering reactions such as phosphates.

KIT COMPONENTS

| Reagents: | REDOXBFP | REDOXBFP20 |
|---|----------|------------|
| d-ROMs Test FAST: | | |
| - Reagent R1 d-ROMs Test FAST Chromogenic mixture condensed in the cuvette, pre-dosed. Store at 15-25°C | 50 pcs | 20 pcs |
| - Reagent R2 d-ROMs Test FAST Buffer pH 4,8, preservatives and stabilizers in ready prepared micro test tubes. Store at 15-25°C | 50 pcs | 20 pcs |
| - Reagent R3 d-ROMs Test FAST Catalyst solution Store at 15-25 °C | 1x1,5ml | 1x0,8ml |
| PAT TEST: | 50 pcs | 20 pcs |
| - Reagent R1 PAT TEST Chromogenic mixture in the cuvette, pre-dosed. Store at 15-25°C | 1x3ml | 1x1,2 ml |
| - Reagent R2 PAT TEST Ferric nitrate solution, stabilizers and preservatives. Store at 15-25°C | | |
| Materials: | 200 pcs | 80 pcs |
| - tips, disposable | 50 pcs | 20 pcs |
| - sterile lancets, disposable | 50 pcs | 20 pcs |
| - heparinized microvettes, disposable | | |

MATERIAL REQUIRED BUT NOT SUPPLIED

Pipette (10 µl and 40 µl) - Photometer

PRECAUTIONS AND WARNINGS

According to Reg 1272/2008 CLP the product is not classified as dangerous.

Handle the product with caution, according to good laboratory practice, avoid ingestion, avoid contact with skin, eyes and mucous membranes. Safety data sheets of individual components are available upon request.

Treat all the samples as if they contain infectious agents. Observe precautionary practices against microbiological hazards by following the procedures and standards for proper sample disposal.

Do not use after the expiration date.

SPECIMENBiological fluid: heparinized capillary blood from which we obtain plasma by centrifugation. **Do not use blood treated with citrate or EDTA** or other iron chelators. Only use ethyl alcohol as a disinfectant.

Storage of the sample:

| 18-25°C | 0-4°C | -20°C |
|----------|----------|----------|
| 12 hours | 48 hours | 2 months |

The test must be performed on the patient fasting from the previous evening. The food and drinks going to change the level of antioxidant in systematically and therefore distort the test results. **Do not** occur substantial change in the value over if you do not vary the condition of the subject.**SAMPLE PREPARATION**

Prepare a microvette by pulling it out from its container and pulling the small cap attached to the main cap.

- Gently massage the fingertip and disinfect it with alcohol (strictly avoid hydrogen peroxide, a powerful oxidant).

- Practice a puncture on a fingertip with a sterile lancet and gently massage the finger in order to promote blood spill.

- Wipe away the first drop of blood (rich in cellular liquid) with a cotton ball.

- Move the finger close to the microvette and let the blood enter through the smaller hole. Fill the microvette until the end of the fins.

- Close the small hole with the cap FIRST and THEN the main cap. Finally insert the microvette into the container.

TEST PROCEDURE**ATTENTION:** Before carrying out the tests, you should prepare all necessary materials. The instrument must be switched on at least 10 minutes before a test and should not be used until it has finished warming up.**d-ROMs test FAST:**

- Centrifuge the capillary blood taken, in order to separate the plasma from the rest of the sample. Always remember to balance the sample during centrifugation.

- Preparer working solution depositing in micro test tube containing reagent R2 10 ul of reagent R3 using the white pipette and shake for reverse for about 10 seconds.

- After centrifugation and preparation of the working solution taken with the pipette (White) that comes with the tool, 10 ul of plasma and deposited in micro test tubes containing the working solution and mix by inversion for at least 10 seconds.

- Open the cap of the vial and remove the safety ring from the cap by pulling it down. Transfer the contents into a cuvette containing the reagent R1, the chromogen condensed.

- Close the vial with the cap and mix gently by inversion for 10 seconds. Avoid foaming.

- Insert the cuvette into the reading room making sure that the sides are ribbed oriented as indicated on the label

PAT test:

- Take the cuvette containing the reagent R1 and add 40 µl of reagent R2 by using the green pipette supplied with the instrument and its disposable tips.

- Close the cuvette with the cap and mix by inversion for exactly 10 seconds.

- Insert the cuvette in the reading cell of the instrument making sure that the knurled sides are oriented as indicated by the label on the instrument. The instrument completes the first reading in about 2 seconds. At the end remove the cuvette from the reading cell.

REFERENCE VALUES**d-ROMs Test FAST:**

| | |
|-----------------------|----------------------------|
| >500 U Carr | Very high oxidative stress |
| 400-500 U Carr | High oxidative stress |
| 340-400 U Carr | Oxidative stress |
| 320-340 U Carr | Slight oxidative stress |
| 300-320 U Carr | Borderline range |
| 250-300 U Carr | Normal values |

PAT TEST:

| | |
|------------------------|--------------------------|
| >2800 U Cor | Very high value |
| 2800-2200 U Cor | Normal value |
| 2200-2000 U Cor | Borderline range |
| 2000-1800 U Cor | Slight deficiency status |
| <1800 U Cor | Deficiency status |

ANALYTICAL PERFORMANCE**d-ROMs Test FAST:**

Unit of measure: 1 U Carr = 0.08 mg / dl hydrogen peroxide

Linearity: the method is linear in the 50-600 U Carr range

Precision: CV% <5.0%

Interference: the addition of anticoagulants capable of chelating iron, such as EDTA or citrates, give rise to underestimates of the data; the use of disinfectants different by ethyl alcohol can lead to anomalous results.

PAT TEST:

Unit of measure: 1 U Cor = 1.4 µmol / liter of vitamin C

Linearity: the method is linear in the 500-10,000 U Cor range

Precision: CV% <5.5%

Interferences: no interference was observed in the presence of a phosphate concentration of less than 40 mg/dl; the addition of anticoagulants capable of chelating iron, such as EDTA or citrates, give rise to underestimates of the data; the use of disinfectants different by ethyl alcohol can lead to anomalous results.

WASTE DISPOSAL

The product should be disposed of according to local legislation on waste management

BIBLIOGRAPHYCornelli U, et al. *Intern. Union of Angiology's Bulletin*. 1999. 15: 7-10.Cesarone MR, et al. *International Angiology*. 1999. 18 (2): 127-130.Alberti A, et al. *Res Chem Intermed*. 2000. 26 (3): 253-67.Trotti R, et al. 2001. *Haematologica*. 86: 85-91.Gerardi GM, et al. *Clin Chem Lab Med*. 2002. 40 (2): 104-110.**SYMBOLS**

| | | | |
|--|---|------------|----------------------------------|
| | Check the instructions for use | REF | Catalog number |
| | Temperature range in which to store the product | LOT | Lot number |
| | Keep away from direct light sources | IVD | For in vitro diagnostic use only |
| | Manufacturer | | Expiry date (year / month) |
| | Disposable | | |

Please read this leaflet carefully before proceeding to use the product. No liability is accepted for damage resulting from improper use or use not covered in this package insert.

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